

PLD5 Antibody



PACO38554

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

PLD5

Uniprot

Q8N7P1

Synonyms:

Inactive phospholipase D5 (Inactive PLD 5) (Inactive choline phosphatase 5) (Inactive phosphatidylcholine-hydrolyzing phospholipase D5) (PLDc), PLD5


Immunogen:

Recombinant Human Inactive phospholipase D5 protein (92-536AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO38554 at dilution of 1:100.