

SETD9 Antibody



PACO38578

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:**Gene ID:**

SETD9

Uniprot

Q8NE22

Synonyms:

SET domain-containing protein 9 (EC 2.1.1. -), SETD9, C5orf35


Immunogen:

Recombinant Human SET domain-containing protein 9 protein (1-299AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: SETD9 antibody at 6 μ g/ml + K562 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 35, 31 kDa. Observed band size: 35 kDa.