

C15orf40 Antibody



PACO38622

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C15orf40

Uniprot

Q8WUR7

Synonyms:

UPF0235 protein C15orf40, C15orf40


Immunogen:

Recombinant Human UPF0235 protein C15orf40 protein (28-153AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO38622 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skin tissue using PACO38622 at dilution of 1:100.