

C15orf40 Antibody, FITC conjugated



PACO38624

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	C15orf40
Source:	Uniprot
Rabbit	Q8WUR7
Isotype:	Synonyms:
IgG	UPF0235 protein C15orf40, C15orf40
Applications:	Immunogen:
ELISA	Recombinant Human UPF0235 protein C15orf40 protein (28-153AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A