

C10orf82 Antibody, FITC conjugated



PACO38640

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

C10orf82

Uniprot

Q8WW14

Synonyms:

Uncharacterized protein C10orf82, C10orf82

Immunogen:

Recombinant Human Uncharacterized protein C10orf82 protein (1-154AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A