

N4BP2L2 Antibody



PACO38662

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

extracellular exosome, nucleus, transcriptional repressor complex, enzyme binding, RNA polymerase II transcription corepressor activity, negative regulation of hematopoietic stem cell differentiation, negative regulation of transcription from RNA polymerase II promoter, positive regulation of hematopoietic stem cell proliferation.

Gene ID:

N4BP2L2

Uniprot

Q92802

Synonyms:

NEDD4-binding protein 2-like 2 (Phosphonoformate immuno-associated protein 5), N4BP2L2, CG005 PFAAP5


Immunogen:

Recombinant Human NEDD4-binding protein 2-like 2 protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO38662 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO38662 at dilution of 1:100.