

PACO38698

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytosol, uridine kinase activity, pyrimidine nucleobase metabolic process, pyrimidine nucleoside salvage.

Gene ID:

UPRT

Uniprot

Q96BW1

Synonyms:

Uracil phosphoribosyltransferase homolog, UPRT


Immunogen:

Recombinant Human Uracil phosphoribosyltransferase homolog protein (1-309AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO38698 at dilution of 1:100.