

# VSIG2 Antibody



PACO38734

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

VSIG2

**Uniprot**

Q96IQ7

**Synonyms:**

V-set and immunoglobulin domain-containing protein 2 (Cortical thymocyte-like protein) (CT-like protein), VSIG2, CTH CTXL

**Immunogen:**


Recombinant Human V-set and immunoglobulin domain-containing protein 2 protein (24-243AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human breast cancer using PACO38734 at dilution of 1:100.