

# STAC3 Antibody



PACO38754

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000

**Protein Background:**

In skeletal muscles contraction, may play a role in neuromuscular synaptic transmission.

**Gene ID:**

STAC3

**Uniprot**

Q96MF2

**Synonyms:**

SH3 and cysteine-rich domain-containing protein 3, STAC3

**Immunogen:**


Recombinant Human SH3 and cysteine-rich domain-containing protein 3 protein (1-253AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western blot. All lanes: STAC3 antibody at 12 $\mu$ g/ml + Rat lung tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 42, 38, 21 kDa. Observed band size: 42 kDa.