

STAC3 Antibody, HRP conjugated



PACO38755

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

In skeletal muscles contraction, may play a role in neuromuscular synaptic transmission.

Gene ID:

STAC3

Uniprot

Q96MF2

Synonyms:

SH3 and cysteine-rich domain-containing protein 3, STAC3

Immunogen:

Recombinant Human SH3 and cysteine-rich domain-containing protein 3 protein (1-253AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A