

PACO38766

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Probable transcriptional regulator involved in the development of the skeleton of the limb, vertebral column and head. Acts by controlling the number of mesenchymal precursor cells and chondrocytes.

Gene ID:

TBX15

Uniprot

Q96SF7

Synonyms:

T-box transcription factor TBX15 (T-box protein 15) (T-box transcription factor TBX14) (T-box protein 14), TBX15, TBX14


Immunogen:

Recombinant Human T-box transcription factor TBX15 protein (1-496AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: TBX15 antibody at 5 μ g/ml + HeLa whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 66, 55 kDa. Observed band size: 66 kDa.