

# SRFBP1 Antibody, Biotin conjugated



PACO38821

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in regulating transcriptional activation of cardiac genes during the aging process. May play a role in biosynthesis and/or processing of SLC2A4 in adipose cells.

**Gene ID:**

SRFBP1

**Uniprot**

Q8NEF9

**Synonyms:**

Serum response factor-binding protein 1 (SRF-dependent transcription regulation-associated protein) (p49/STRAP), SRFBP1

**Immunogen:**

Recombinant Human Serum response factor-binding protein 1 protein (201-429AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A