

RFESD Antibody



PACO38826

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

RFESD

Uniprot

Q8TAC1

Synonyms:

Rieske domain-containing protein, RFESD


Immunogen:

Recombinant Human Rieske domain-containing protein (1-157AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO38826 at dilution of 1:100.