

# SLC25A40 Antibody, HRP conjugated



PACO38831

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SLC25A40

**Uniprot**

Q8TBP6

**Synonyms:**

Solute carrier family 25 member 40 (Mitochondrial carrier family protein), SLC25A40, MCFP

**Immunogen:**

Recombinant Human Solute carrier family 25 member 40 protein (40-102AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A