

DEFB104A Antibody



PACO38846

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Has antimicrobial activity. Synergistic effects with lysozyme and DEFB103.

Gene ID:

DEFB104A

Uniprot

Q8WTQ1

Synonyms:

Beta-defensin 104 (Beta-defensin 4) (BD-4) (DEFB-4) (hBD-4) (Defensin, beta 104), DEFB104A; DEFB104B, DEFB104 DEFB4;


Immunogen:

Recombinant Human Beta-defensin 104 protein (23-72AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human melanoma using PACO38846 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO38846 at dilution of 1:100.