

CRYGN Antibody



PACO38878

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

CRYGN

Uniprot

Q8WXF5

Synonyms:

Gamma-crystallin N (Gamma-N-crystallin), CRYGN


Immunogen:

Recombinant Human Gamma-crystallin N protein (1-182AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO38878 at dilution of 1:100.