

PACO38910

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

nuclear chromatin, chromatin DNA binding, RNA polymerase II core promoter sequence-specific DNA binding, RNA polymerase II transcription coactivator activity, cellular response to estrogen stimulus, establishment of protein localization, positive regulation of gene expression, epigenetic, positive regulation of histone H3-K14 acetylation, positive regulation of transcription from RNA polymerase II promoter.

**Gene ID:**

WBP2

**Uniprot**

Q969T9

**Synonyms:**

WW domain-binding protein 2 (WBP-2), WBP2

**Immunogen:**


Recombinant Human WW domain-binding protein 2 protein (1-261AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO38910 at dilution of 1:100.