

# RAB37 Antibody



PACO38922

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

endoplasmic reticulum-Golgi intermediate compartment, endosome, plasma membrane, secretory granule membrane, GTPase activity, intracellular protein transport, protein localization to plasma membrane, protein secretion, Rab protein signal transduction, regulation of exocytosis.

**Gene ID:**

RAB37

**Uniprot**

Q96AX2

**Synonyms:**

Ras-related protein Rab-37, RAB37

**Immunogen:**


Recombinant Human Ras-related protein Rab-37 protein (2-220AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human skin tissue using PACO38922 at dilution of 1:100.