

CXorf40B Antibody



PACO38946

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CXorf40B

Uniprot

Q96DE9

Synonyms:

Protein CXorf40B, CXorf40B

Immunogen:

Recombinant Human Protein CXorf40B protein (1-158AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A