

CCDC42 Antibody



PACO39006

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

CCDC42

Uniprot

Q96M95

Synonyms:

Coiled-coil domain-containing protein 42, CCDC42, CCDC42A


Immunogen:

Recombinant Human Coiled-coil domain-containing protein 42A protein (1-242AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: CCDC42 antibody at 0.8 μ g/ml. Lane 1: Mouse brain tissue. Lane 2: Mouse heart tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 39, 30 kDa. Observed band size: 39 kDa.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39006 at dilution of 1:100.