

# TSGA13 Antibody



PACO39038

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000,  
IHC:1:20-1:200

**Protein Background:**

positive regulation of protein targeting to mitochondrion.

**Gene ID:**

TSGA13

**Uniprot**

Q96PP4

**Synonyms:**

Testis-specific gene 13 protein, TSGA13

**Immunogen:**

Recombinant Human Testis-specific gene 13 protein (1-275AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images


---



Western blot. All lanes: TSGA13 antibody at 4 $\mu$ g/ml + A549 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 32 kDa. Observed band size: 32 kDa.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO39038 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39038 at dilution of 1:100.