

LYG1 Antibody



PACO39066

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

LYG1

Uniprot

Q8N1E2

Synonyms:

Lysozyme g-like protein 1 (EC 3.2.1. -), LYG1


Immunogen:

Recombinant Human Lysozyme g-like protein 1 protein (20-194AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO39066 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon tissue using PACO39066 at dilution of 1:100.