

PAF1 Antibody



PACO39102

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:2000-1:10000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

PAF1

Uniprot

Q8N7H5

Synonyms:

RNA polymerase II-associated factor 1 homolog (hPAF1) (Pancreatic differentiation protein 2), PAF1, PD2


Immunogen:

Recombinant Human RNA polymerase II-associated factor 1 homolog protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: HepG2 whole cell lysate, HeLa whole cell lysate, K562 whole cell lysate, Mouse brain tissue. All lanes: PAF1 antibody at 3.2 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 60, 56 kDa. Observed band size: 60, 48 kDa.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39102 at dilution of 1:100.