

C8orf31 Antibody



PACO39110

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

C8orf31

Uniprot

Q8N9H6

Synonyms:

Uncharacterized protein C8orf31, C8orf31

Immunogen:

Recombinant Human Uncharacterized protein C8orf31 protein (1-62AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A