

# FAM71C Antibody



PACO39146

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000

**Protein Background:**

nucleus.

**Gene ID:**

FAM71C

**Uniprot**

Q8NEG0

**Synonyms:**

Protein FAM71C, FAM71C

**Immunogen:**


Recombinant Human Protein FAM71C protein (1-241AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western blot. All lanes: FAM71C antibody at 6 $\mu$ g/ml. Lane 1: 293T whole cell lysate. Lane 2: Jurkat whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 28 kDa. Observed band size: 28 kDa.