

C1orf64 Antibody



PACO39150

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SRARP

Uniprot

Q8NEQ6

Synonyms:

Steroid receptor-associated and regulated protein (ER-related factor) (Steroid receptor-regulated protein), SRARP, C1orf64 ERRF SSPR

Immunogen:

Recombinant Human Uncharacterized protein C1orf64 protein (1-169AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A