

PACO39174

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Hormone that specifically promotes pancreatic beta cell proliferation and beta cell mass expansion, thereby improving glucose tolerance. Promotes pancreatic beta cell proliferation without insulin resistance. Also acts as a blood lipid regulator by regulating serum triglyceride levels; possibly by promoting ANGPTL3 cleavage.

Gene ID:

Angptl8

Uniprot

Q8R1L8

Synonyms:

Angiopoietin-like protein 8 (Betatrophin) (Lipasin) (Refeeding-induced fat and liver protein), Angptl8, Gm6484 R1fl

Immunogen:

Recombinant Mouse Betatrophin protein (16-198AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A