

LMO3 Antibody



PACO39178

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

cytoplasm, negative regulation of ERK1 and ERK2 cascade, positive regulation of fat cell differentiation, positive regulation of glucocorticoid receptor signaling pathway, positive regulation of peroxisome proliferator activated receptor signaling pathway.

Gene ID:

LMO3

Uniprot

Q8TAP4

Synonyms:

LIM domain only protein 3 (LMO-3) (Neuronal-specific transcription factor DAT1) (Rhombotin-3), LMO3, RBTN3 RBTNL2 RHOM3


Immunogen:

Recombinant Human LIM domain only protein 3 protein (1-145AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: LMO3 antibody at 6 μ g/ml + HeLa whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 17, 20, 18, 19 kDa. Observed band size: 17 kDa.



Immunohistochemistry of paraffin-embedded human colon tissue using PACO39178 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO39178 at dilution of 1:100.