

CMTM2 Antibody, FITC conjugated



PACO39184

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

nucleoplasm.

Gene ID:

CMTM2

Uniprot

Q8TAZ6

Synonyms:

CKLF-like MARVEL transmembrane domain-containing protein 2 (Chemokine-like factor superfamily member 2), CMTM2, CKLFSF2

Immunogen:

Recombinant Human CKLF-like MARVEL transmembrane domain-containing protein 2 protein (1-80AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A