

# SSMEM1 Antibody, Biotin conjugated



PACO39221

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

SSMEM1

**Uniprot**

Q8WWF3

**Synonyms:**

Serine-rich single-pass membrane protein 1, SSMEM1, C7orf45

**Immunogen:**

Recombinant Human Serine-rich single-pass membrane protein 1 protein (56-244AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A