

# Scgb3a2 Antibody



PACO39234

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

extracellular region.

**Gene ID:**

Scgb3a2

**Uniprot**

Q920H1

**Synonyms:**

Secretoglobin family 3A member 2 (Pneumo secretory protein 1) (PnSP-1) (Uteroglobin-related protein 1), Scgb3a2, Pnsp1 Ugrp1

**Immunogen:**

Recombinant Mouse Secretoglobin family 3A member 2 protein (22-139AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A