

CHMP6 Antibody, FITC conjugated



PACO39320

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CHMP6

Uniprot

Q96FZ7

Synonyms:

Charged multivesicular body protein 6 (Chromatin-modifying protein 6) (Vacuolar protein sorting-associated protein 20) (Vps20) (hVps20), CHMP6, VPS20

Immunogen:

Recombinant Human Charged multivesicular body protein 6 protein (2-201AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A