

SLC22A18AS Antibody, HRP conjugated



PACO39383

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

SLC22A18AS

Uniprot

Q8N1D0

Synonyms:

Beckwith-Wiedemann syndrome chromosomal region 1 candidate gene B protein (Organic cation transporter-like protein 2 antisense protein) (Solute carrier family 22 member 1-like antisense protein) (Solute carrier family 22 member 18 antisense protein) (p27-Beckwith-Wiedemann region 1 B) (p27-BWR1B), SLC22A18AS, BWR1B, BWSCR1B, ORCTL2S, SLC22A1LS

Immunogen:

Recombinant Human Beckwith-Wiedemann syndrome chromosomal region 1 candidate gene B protein (1-150AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A