

STPG2 Antibody



PACO39402

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

STPG2

Uniprot

Q8N412

Synonyms:

Sperm-tail PG-rich repeat-containing protein 2, STPG2, C4orf37


Immunogen:

Recombinant Human Sperm-tail PG-rich repeat-containing protein 2 protein (201-459AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: STPG2 antibody at 0.8 μ g/ml + HL60 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 51 kDa. Observed band size: 51 kDa.



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO39402 at dilution of 1:100.