

C17orf78 Antibody



PACO39406

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:**Gene ID:**

C17orf78

Uniprot

Q8N4C9

Synonyms:

Uncharacterized protein C17orf78, C17orf78


Immunogen:

Recombinant Human Uncharacterized protein C17orf78 protein (1-186AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO39406 at dilution of 1:100.



Immunofluorescent analysis of A549 cells using PACO39406 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).