

CCDC82 Antibody



PACO39410

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CCDC82

Uniprot

Q8N4S0

Synonyms:

Coiled-coil domain-containing protein 82, CCDC82

Immunogen:

Recombinant Human Coiled-coil domain-containing protein 82 protein (1-344AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A