

FAM71E2 Antibody



PACO39422

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

nucleus.

Gene ID:

FAM71E2

Uniprot

Q8N5Q1

Synonyms:

Protein FAM71E2, FAM71E2, C19orf16

Immunogen:

Recombinant Human Protein FAM71E2 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A