

ZFAND2A Antibody



PACO39430

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

ZFAND2A

Uniprot

Q8N6M9

Synonyms:

AN1-type zinc finger protein 2A, ZFAND2A


Immunogen:

Recombinant Human AN1-type zinc finger protein 2A protein (1-145AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO39430 at dilution of 1:100.