

PACO39486

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Required for calcium-mediated regulation of ciliary beat frequency and mucociliary clearance in the airway. Might be involved in the regulation of intracellular calcium in tracheal epithelial cells.

**Gene ID:**

AGR3

**Uniprot**

Q8TD06

**Synonyms:**

Anterior gradient protein 3 (AG-3) (AG3) (hAG-3) (Anterior gradient 3 homolog) (Breast cancer membrane protein 11) (Protein disulfide isomerase family A, member 18), AGR3, BCMP11 PDIA18

**Immunogen:**

Recombinant Human Anterior gradient protein 3 protein (22-166AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A