

GPX8 Antibody, HRP conjugated



PACO39491

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

endoplasmic reticulum lumen, peroxidase activity, response to reactive oxygen species.

Gene ID:

GPX8

Uniprot

Q8TED1

Synonyms:

Probable glutathione peroxidase 8 (GPx-8) (GSHPx-8) (EC 1.11.1.9), GPX8

Immunogen:

Recombinant Human Probable glutathione peroxidase 8 protein (41-209AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A