

ARL11 Antibody, HRP conjugated



PACO39523

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in apoptosis. May act as a tumor suppressor.

Gene ID:

ARL11

Uniprot

Q969Q4

Synonyms:

ADP-ribosylation factor-like protein 11 (ADP-ribosylation factor-like tumor suppressor protein 1), ARL11, ARLTS1

Immunogen:

Recombinant Human ADP-ribosylation factor-like protein 11 protein (2-196AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A