

PACO39626

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

UBL7

Uniprot

Q96S82

Synonyms:

Ubiquitin-like protein 7 (Bone marrow stromal cell ubiquitin-like protein) (BMSC-UbP)
(Ubiquitin-like protein SB132), UBL7, BMSCUBP


Immunogen:

Recombinant Human Ubiquitin-like protein 7 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO39626 at dilution of 1:100.