

SLC22A24 Antibody, FITC conjugated



PACO39648

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

integral component of plasma membrane, inorganic anion exchanger activity, sodium-independent organic anion transmembrane transporter activity, anion transmembrane transport, sodium-independent organic anion transport.

Gene ID:

SLC22A24

Uniprot

Q8N4F4

Synonyms:

Solute carrier family 22 member 24, SLC22A24

Immunogen:

Recombinant Human Solute carrier family 22 member 24 protein (37-145AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A