

CD177 Antibody, HRP conjugated



PACO39655

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular exosome, plasma membrane, blood coagulation, leukocyte migration.

Gene ID:

CD177

Uniprot

Q8N6Q3

Synonyms:

CD177 antigen (Human neutrophil alloantigen 2a) (HNA-2a) (NB1 glycoprotein) (NB1 GP) (Polycythemia rubra vera protein 1) (PRV-1) (CD antigen CD177), CD177, NB1 PRV1

Immunogen:

Recombinant Human CD177 antigen protein (22-321AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A