

# MLKL Antibody, FITC conjugated



PACO39684

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Pseudokinase that plays a key role in TNF-induced necroptosis, a programmed cell death process. Activated following phosphorylation by RIPK3, leading to homotrimerization, localization to the plasma membrane and execution of programmed necrosis characterized by calcium influx and plasma membrane damage. Does not have protein kinase activity.

**Gene ID:**

MLKL

**Uniprot**

Q8NB16

**Synonyms:**

Mixed lineage kinase domain-like protein (hMLKL), MLKL

**Immunogen:**

Recombinant Human Mixed lineage kinase domain-like protein (1-471AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A