

MLKL Antibody, Biotin conjugated



PACO39685

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Pseudokinase that plays a key role in TNF-induced necroptosis, a programmed cell death process. Activated following phosphorylation by RIPK3, leading to homotrimerization, localization to the plasma membrane and execution of programmed necrosis characterized by calcium influx and plasma membrane damage. Does not have protein kinase activity.

Gene ID:

MLKL

Uniprot

Q8NB16

Synonyms:

Mixed lineage kinase domain-like protein (hMLKL), MLKL

Immunogen:

Recombinant Human Mixed lineage kinase domain-like protein (1-471AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A