

C17orf47 Antibody




PACO39698


Product Information

Size: 50ug	Protein Background:
Reactivity: Human	Gene ID: C17orf47
Source: Rabbit	Uniprot Q8NEP4
Isotype: IgG	Synonyms: Uncharacterized protein C17orf47, C17orf47
Applications: ELISA, IHC	Immunogen: Recombinant Human Uncharacterized protein C17orf47 protein (1-300AA).
Recommended dilutions: ELISA:1:2000-1:10000, IHC:1:20-1:200	Storage: Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO39698 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39698 at dilution of 1:100.