

NFAM1 Antibody, Biotin conjugated



PACO39705

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May function in immune system as a receptor which activates via the calcineurin/NFAT-signaling pathway the downstream cytokine gene promoters. Activates the transcription of IL-13 and TNF-alpha promoters. May be involved in the regulation of B-cell, but not T-cell, development. Overexpression activates downstream effectors without ligand binding or antibody cross-linking.

Gene ID:

NFAM1

Uniprot

Q8NET5

Synonyms:

NFAT activation molecule 1 (Calcineurin/NFAT-activating ITAM-containing protein) (NFAT-activating protein with ITAM motif 1), NFAM1, CNAIP

Immunogen:

Recombinant Human NFAT activation molecule 1 protein (186-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A