

FCGRT Antibody



PACO39714

Product Information

Size:

50ug

Reactivity:

Macaca fascicularis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds to the Fc region of monomeric immunoglobulins gamma. Mediates the uptake of IgG from milk. Possible role in transfer of immunoglobulin G from mother to fetus.

Gene ID:

FCGRT

Uniprot

Q8SPV9

Synonyms:

IgG receptor FcRn large subunit p51 (FcRn) (IgG Fc fragment receptor transporter alpha chain) (Neonatal Fc receptor), FCGRT, FCRN

Immunogen:

Recombinant Macaca fascicularis IgG receptor FcRn large subunit p51 protein (24-297AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A