

TRUB1 Antibody



PACO39766

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be responsible for synthesis of pseudouridine from uracil in transfer RNAs.

Gene ID:

TRUB1

Uniprot

Q8WWH5

Synonyms:

Probable tRNA pseudouridine synthase 1 (EC 5.4.99. -), TRUB1, PUS4

Immunogen:

Recombinant Human Probable tRNA pseudouridine synthase 1 protein (2-349AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A